**南通海发制药设备有限公司**  
High Fine Pharmaceutical Equipment Co., Ltd  
**南通海发水处理工程有限公司**  
High Fine Engineering Ltd

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*FDA standard water system  
Mainstream supplier of global water system*

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Water system equipment that meet FDA standards

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Main supplier of global water system equipment

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
HIGHFINE always takes the economic benefits of pharmaceutical companies, social responsibilities and human health as our value of survival. We would like to do diligent work, to do what we can and continuous dedication of life in order to realize health undertakings that can invigorate a nation with technology and that can bring benefit to the people in Chinese hope land together with the friends who support and care about the industrialization of China pharmaceutical industry.

携手共进 **26**年  
 26 YEARS, HAND IN HAND  
 COMMON DEVELOPMENT

A collection of logos for various pharmaceutical and chemical companies, arranged in a grid-like pattern across the bottom half of the page. The logos include:

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HIGHFINE mainly provides the following equipments for domestic and foreign main pharmaceutical enterprises:

Design, installation and verification of pharmaceutical water preparation units  
制药用水制备机组的设计、安装和验证

RO EDI纯化水制备机组  
RO EDI purified water preparation unit

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Design, installation and verification of distribution unit  
分配机组的设计、安装和验证

PW纯化水分配机组  
PW purified water distribution unit

WFI注射用水分配机组  
WFI water for injection distribution unit

Design, manufacturing and installation of mixing liquid tank and liquid conveying system  
配料罐及配料系统设计、制造和安装

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High Fine is the only joint venture in China to manufacture pharmaceutical water equipments professionally. The Company adopts the European and American technologies to introduce the world leading manufacturing technology into China and to provide water system equipment that completely meet the requirements of European and American pharmacopoeia. In addition, the Company has international authoritative engineers in pharmaceutical water equipment field. Under this circumstance, the design of each set of equipment, the equipment manufacturing, construction, the quality of authentication service and after-sale service are assured.

HIGHFINE is not only the mainstream supplier of the water system equipments in Chinese pharmaceutical industry. Our struggle objective and the ultimate ideal are to become the leader and mainstream supplier of the water system equipments in global pharmaceutical industry.

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
Time passed quickly. We from young to mature, from full of enthusiasm in the past into today's more persistent and more rigorous, because we deeply feel that the responsibility be weightier than Mount Tai. Time made our ideal specific, given its true meaning. If human's living environment better improved as our efforts, then we will feel that HighFine's wisdom is infinite.

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Engineering design consulting and management verification

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Simply speaking, the responsibility of validation is assumed by developer (design), quality assurance department (certify) and users (operate). The responsibility of developer is providing system development and technology support, compiling and making necessary test for documents. While the quality assurance department is responsible for constituting validation guideline and steps.

A process which is correct and had been verified can sequentially produce products which meet with intending cost and quality and are of high purity and safety. Then they can reduce cost, guarantee quality and comply with GMP specification and FDA requirement.

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Our company provide following validation consulting service:

- establishing validation master plan
- compiling validation scheme
- \*Design qualification (DQ)      \*Installation qualification (IQ)
- \*Calibration      \*Operation qualification (OQ)
- \*Performance qualification (PQ)
- Computer system validation
- Compiling standard operation procedure (SOP)
- Training

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## Fine Management

Every tiny detail is worthy of our close attention, every small action is the part which HighFine's staff working in. Only to give enough attention, correct tiny deviation, can the enterprise achieve success. HighFine based on this standard, and constantly strive to improve ourselves to be perfect as our duty.

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Automatic longitudinal seam welding machine

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Automatic circumferential seam welding machine

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Tube plate automatic welding machine

Orbitalum GF4切管机

ORBITALUM GF4 pipe cutting machine

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Automatic track welding

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Laser cutting machine

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Gantry CNC boring and milling machine

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## Lean production

HighFine's lean production management system includes 4 parts. They are the cost quota control system, the work hours quota management system, process menu system, 6s on-site management system.

Through strict control in the production process of consumption of materials, energy, artificial time, process design, site management requirements, to realize the target: the lean production management system with minimization of cost and benefit maximization.

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Specialization, Securitization, Standardization

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
## The unique "3S" quality control system

HighFine QC testing center.  
To provide efficient reliable testing and quality control for each product.  
To ensure that each procedure to meet the requirements of customers.  
We do our best to make customers' satisfaction.

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► Purified water system

We take modular design as the concept for the purified water system design to make the purified water system has such characteristics of independence, modularization and simple operation, etc, which endures the consecutiveness, stability and reliability of equipment operation. In addition, we will make more than two design plans for the selection by the clients in accordance with the water quality report and URS provided by the clients, including pretreatment configuration, combination of desalination system, selection of distribution system and disinfection methods, etc. In the meantime, we also provide a complete set of verification system for clients, including design validation in the preliminary stage, FAT test within the factory, installation validation in later stage, operation validation and performance validation as well as perfect after-service system.

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► Clean steam generator


The clean steam usually takes purified water as feeding water and the steam is produced through pure steam generator or the first effect evaporator of multi-effect water distiller. At the time of condensation, the clean steam shall meet the requirements of water for injection. Deionized water and purified water can all be taken as feeding water of clean steam generator. After evaporation and separation (remove such pollutants as particles and bacterial endotoxin), it can be conveyed to the using point under certain pressure.

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► Multi-effect water distiller

Distillation is the technical process for chemical and microbiological purification of material water through steam-liquid phase transition method and separation method. During this technical process, the water is evaporated and the steam produced isolated from the water. What's more, the unevaporated water that flow back dissolved solid substance, non-volatile matter and high molecular weight impurities. During the process of distillation, low molecular impurities may be taken along in the steam produced by water evaporation in the form of water spray or water drop. Under this circumstance, we need a separating unit to remove the thin water spray and the carried impurities, including endotoxin. The purified steam after condensation will become water for injection. At least 99.99% endotoxin content may decrease through distillation method.

Systematic · overall · professional  
product classification

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### ► Scope of application:

Pure water in such industries as electronics, medicine and food, etc; purification and preparation for water used in textile industry and chemical industry; purification and preparation for water used in food and drinks and enological process; concentration and recovery of the useful materials of aqueous solution in industrial production; preliminary desalination processing of high pressure boiler make-up water of the enterprises, such as electrical power plant, etc.; brackish water and seawater desalination; as the first class desalinating equipment for high purity water production.

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Water system equipment / RO + EDI system

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EDI equipment belongs to fine treatment water system and usually used in cooperation with reverse osmosis (RO) to constitute the ultrapure water treatment system with pretreatment, reverse osmosis and EDI devices. EDI equipment replaced the mixed ion-exchange unit used in traditional water treatment technology. EDI equipment requires the electrical resistivity of the water to reach 0.025-0.5MΩ·cm and the reverse osmosis device can completely meet the requirements. EDI equipment can produce ultrapure water with the electrical resistivity above 15MΩ·cm.

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Water system equipment  
Hot water disinfection purification water system

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### ▶ Pure water equipment requirements:

- The water quality meets all the provisions of China pharmacopoeia standard and FDA;
- Automatic equipment operation and conditional automatic processing program (such as back flush, regeneration, acid-washing and disinfection procedures);
- Monomer and pipe equipment meet GMP requirements (all such back end processing equipments as sterilizer, membrane filtration, terminal water tank and pipeline adopt the 316L material and the pre-treatment pipeline adopts SUS304 material).

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Water system equipment / Second class RO+EDI purified water system

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- 电导率 $\leq 0.5\mu\text{s}/\text{cm}$
- $\text{TOC} \leq 200\text{ppb}$
- 微生物 $\leq 20\text{cfu}/\text{ml}$

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Multi-effect water distiller  
S series multi-effect water distiller

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S series multi-effect water distiller is designed according to the principle of falling film and multiple-stage evaporation. The first effect is to conduct heating by external heat source (industrial steam) and then each pillar can get enough heat energy by using the pure vapor produced from the fore-effect. The material water that did not evaporated in the first effect can be taken as the material in the next effect and so on.

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### ▶ Equipment features:

- The whole equipment is made by 316L stainless steel and the operation surface is conducted with acid-washing and passivating treatment;
- The pipeline system adopts the mirror polish 0.6μmRa 316L stainless steel tube and conducts with automatic track welding and connects with clean flange and pharmaceutical level gasket.
- Use material wool to keep warm and the exterior of the equipment is covered with 314 stainless steel matte finishing protective layer and its scaffold is made of stainless steel.

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Multi-effect water distiller  
Technical parameters of S series multi-effect water distiller

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Inlet pressure of industrial steam is 0.3MPa and the inflow working pressure of cooling water is 0.2MPa

规格型号 Product model	效数 Number of effect	生产能力 Lt/h Production capacity	原料水 Kg/h Material water	工业蒸汽 Kg/h Industrial steam	冷却水 Lt/h Cooling water	高 mm Height	长 mm Length	宽 mm Width
MS100	4	100	115	45	80	2650	1300	710
MS200	4	200	230	76	120	2650	1360	710
MS200	5	200	230	62	50	2650	1620	710
MS300	4	300	345	109	190	2650	1410	730
MS300	5	300	345	89	90	2650	1730	730
MS500	4	500	575	185	310	3170	1680	900
MS500	5	500	575	148	130	3170	2000	900
MS750	4	750	865	256	300	3260	1680	900
MS750	5	750	865	205	178	3260	2000	900
MS1000	5	1000	1150	295	256	3685	2180	980
MS1000	6	1000	1150	245	65	3685	2830	980
MS1500	5	1500	1725	460	390	3685	2680	1100
MS1500	6	1500	1725	380	95	3685	3150	1100
MS2000	5	2000	2300	540	345	3940	3150	1100
MS2000	6	2000	2300	490	130	3950	3580	1100
MS3000	6	3000	3450	730	190	4050	4000	1200
MS4000	6	4000	3450	910	310	4800	4000	1200
MS5000	6	5000	5750	1140	440	4200	4550	1350

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Multi-effect water distiller  
F series multi-effect water distiller

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F series multi-effect water distiller adopts Level 3 centrifugal separation technology. Siemens S7 series control system has Chinese and English operating language, water inlet booster pump with standard configuration, dashpot and collection tank. The equipment is equipped standard with plant steam reducing valve, the device for continuous emission of non-condensable gas. In addition, the F series equipment continuously discharge the residual liquid and record the electrical conductivity of water for injection and the valve location. It is conducted with complete equipment manufacture factory acceptance test and the inlet water uses deionized water or RO water.

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Multi-effect water distiller  
Technical parameters of F series multi-effect water distiller

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Multi-effect water distiller  
B series multi-effect water distiller

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B series multi-effect water distiller adopts rising film and multi-diversion form to evaporate. The industrial steam enters the evaporator shell pass through valves and interchange heat with the purified water rising film in tube pass. The purified water after heat exchange becomes vapour, then impure heat source will be removed from the vapour through diversion separation device in the machine to produce the secondary steam. The pure steam after steam-water separation will become distilled water after being condensed as heating source in II effect. Then the distilled water will enter into the condenser under the pressure differential action and then enter the storage tank through the pipe for storage. The concentrated water without evaporation and separation shall enter the next effect for re-evaporation through concentrated water inlet.

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### ▶ Equipment features:

- Adopts rising film for evaporation to solve the problem of uneven distribution of water.
- High thermal efficiency, 100% heat exchange.
- Save industrial steam consumption, 20% energy saving.
- Distribute feeding water consumption amount for each effect accurately to reduce material water loss.
- No need to exchange the multi-diversion separation device and life-span of the equipment can be extended.

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Multi-effect water distiller  
Technical parameters of B series multi-effect water distiller

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### ▶ Production capacity

The standard equipment model table shows the equipment production capacity under the condition of 3bar pure steam pressure and 6bar industrial steam with the unit in kg/h. The actual production capacity of some type equipments shall be acquired in accordance with the actual operation pressure and use this actual operation pressure to multiply the conversion coefficient in the specified production capacity.

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### ▶ Model selection example

- Industrial steam pressure: 7bar
- Pure steam pressure needed: 2bar
- The highest demand of pure steam: 500kg/h
- Conversion coefficient: 1.7
- Specified output needs:  $500/1.7=294$ kg/h
- Model selected: PSG 300 DTS

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### ▶ Public facility demand

- Industrial saturated steam without impurities and corrosive materials, the maximum pressure is 9bar.
- Feeding water requires deionized water without silica, chlorine, amine and volatile substance, etc. The pressure required is 1-2bar.
- Power supply standards: according to the clients' requirements.
- Dry, oil-free and clean compressed air, the minimum pressure required is 6bar.

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Inlet pressure of industrial steam is 0.3MPa and the pure steam working pressure is 0.2MPa

规格型号 Product model	生产能力 Lt/h Production capacity	原料水 Kg/h Material water	工业蒸汽 Kg/h Industrial steam	高 mm Height	长 mm Length	宽 mm Width
PSG100	100	115	130	2100	900	570
PSG200	200	230	250	2500	900	570
PSG300	300	345	370	2500	1110	680
PSG400	400	460	480	2600	1110	680
PSG500	500	575	173	2600	1170	720
PSG800	800	920	960	3500	1775	930
PSG1000	1000	1150	1200	3700	1775	930
PSG2000	2000	2300	2400	4000	2030	1900
PSG3000	3000	3450	3600	4850	2030	1900
PSG4000	4000	4600	4800	5100	2200	2100

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Storage distribution system  
Pure water storage distribution system

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Distribution system is specially designed according to customers' real situation on the basis of GMP, ISPE and FDA guidelines. The system adopts the frame structure and has such advantages as space-saving, plug-and-play and convenient installation. It will be conducted with pre-installation and operational qualification at the manufacturing plant. All the pipelines use the 316L stainless steel without dead leg and the sanitary type chuck connection mode. The roughness of the surface contacted by all products shall be less than  $0.6\mu\text{m}$ . There is no buffering tank in the loop to avoid bacteria growth, which is in compliance with FDA requirements.

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Storage distribution system  
Distilled water storage distribution system

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Distribution system is specially designed according to customers' real situation on the basis of GMP, ISPE and FDA guidelines. All the pipelines use the 316L stainless steel without dead leg and sanitary type chuck connection mode. The roughness of the surface contacted by all products shall be less than  $0.6\mu\text{m}$ . All the pipe fittings choose ASME-BPE2000 type accessories, which can realize thermostatic control automatically and operate with the water distiller in chain. All the control valves are imported. The whole system can be conducted with the pre-installation and operational qualification at the manufacturing plant and the comprehensive validation documents are provided.

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On-line cleaning system  
Stationary CIP workstation

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CIP cleaning device can be used as the indissoluble production equipment and can also be used in the cleaning system with manual operation and automatic operation. The device is nearly used in such plants as food, drinks and pharmacy, etc. CIP device can not only clean the machine, but can also control the cleaning methods of microorganism.

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On-line cleaning system  
Mobile CIP workstation

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CIP cleaning device has the following advantages:

- Can make the production plan reasonable and improve the production capacity.
- Compared with hand cleaning, CIP does not affect the cleaning effect due to different operators, but can also improve product quality.
- Can avoid the danger in the cleaning operation and save labor force.
- Can save the cost of cleaning agent, steam, water and production.
- Can increase the service life of machine components.
- CIP cleaning devices include such three types as manual operation, semi-automatic and automatic for the selection by the users.

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### Biological products tank equipment in batching system

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Automatic welding equipment that can manufacture stainless steel pressure vessel.

Three kinds of technologies for surface treatment of the stainless steel container:

- 400#+ electrobrightening RaO. 4 $\mu$ m; 400#+600#+800#+2 electrobrightening RaO. 28 $\mu$ m
- Surface sandblasting and surface matte sanding Ra $\leq$ 0.6
- Electrolytic polishing technologies and devices

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Workshop and equipment for mechanical mixing.

Possess the devices for physical and chemical analysis of materials, thickness inspection and inner surface roughness test. The equipment design and manufacturing are totally in compliance with FDA requirements; the pressure vessel meets the national relevant standards and the requirements for pressure vessel safety inspection; meet the relevant requirements for equipment installation.

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HF series bioreactor is the ventilating machine stirred tank reactor and adopts PLC + color touch screen system intelligent control, which can meet the requirements for usage from the lab to pilot plant test and scale of production. This system has such advantages as simple structure, convenient operation and high degree of automation. It is widely used in such fields as edible mushroom, medicinal fungus liquid seeding, biological pharmacy, food additives, organism oil-recovery, biological energy source, biopesticide, bio-fertilizer, biological feed enzyme preparation and vitamin, etc.

- The control parameters mainly include: temperature, rotating speed, ventilatory capacity, tank pressure, pH, Do, feed supplement, defoaming, biomass liveweight and scaling, etc.
- Volume: 5L ~500m<sup>3</sup>

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Our products obtain approval of professional field.

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- Have the design, manufacturing and engineering experience of global pharmaceutical water system
- Have the world first class engineering technology team
- Have the experience for providing services for more than 100 famous clients at home and aboard
- Completely familiar with and grasp the international "Good Manufacturing Practice" and industry guidelines and pharmaceutical industry standard
- Provide the documentary system that can meet FDA, cGMP requirements
- Provide FAT testing documents for complete machine
- Communication and cooperation with the international mainstream research field / grasp the advanced internal technological quality to synchronize with the international